


## **TECHBRIEF**

# **Critical Cleaning for 3D Printing in the Life Sciences**

# Abstract

The technique of three-dimensional (3D) printing can be utilized to create not only simple parts but those with complex geometries as well. These 3D objects are generated via the deposition of successive layers of materials comprising various polymeric and metallic substrates. 3D printing, also known as additive manufacturing, can be easily automated to produce product in commercial quantities. As 3D printing technology has advanced, its use has rapidly increased across a wide range of industries, including life science, healthcare, and medical products, among others.

This TechBrief focuses on the critical cleaning of a wide range of hard surfaces utilized in 3D printing processes across industries, such as extruded plastics, photosensitive substrates, metals, and metal alloys. We use the term critical cleaning to denote situations where the level of cleaning directly impacts the value of the end product or materially increases manufacturing efficiency.



## Wide Variety of Materials

Different resins/materials are preferred for each 3D printing application depending on the 3D printing process and printer that is used, and the design and function of the part being produced. The application may require the material to be clear or colored, rigid or flexible, and/or tough and durable or soft. The material may also need to qualify for use in medical and dental applications. Consequently, biocompatibility, toxicity, biodegradability and other factors of the resins must be considered. It is possible to choose pure or composite materials with specific chemical, fatigue, impact, heat, UV, and water resistance characteristics. There are soluble, elastic, flexible, rigid, and soft materials as well.

## Material Selection for Specific Applications and Printing Approaches

Additive manufacturing is a general term referring to a number of different 3D printing processes that in 2015 were classified into seven categories within the [ISO/ASTM 52900](#) standard: extrusion, vat polymerization, powder bed fusion, material jetting, binder jetting, direct energy deposition, and sheet lamination. The seven categories themselves each encompass many different 3D printing technologies.

The choice of 3D printing process will be dictated by the printer cost and capabilities, the speed of the process, the quality of the parts produced, and the practical implementation of the method. Most importantly, the properties of the product must meet the expectations of the end user.

The choice of material will depend on the type of printing process, the capabilities of the printer, and the nature of the part being fabricated, e.g., complexity, functionality, color, texture, appearance, and so on. In addition, certain filaments are more challenging to work with than others and will likely require more experimentation. The quality of the filament is also important, as chemical contaminants in lower-grade material can have the propensity to clog printer nozzles and lead to poorer-quality parts.

When choosing 3D-printing materials, the first step is to define the performance requirements for the printed part. These requirements can then be translated into appropriate material properties. Examples for plastic filaments include tensile strength, flexural modulus (the tendency of a material to bend), elongation, impact strength, hardness, creep, and compression set. Often, the combined set of material properties will point to just one or two optimal choices.

## Healthcare and Pharmaceutical Applications

### Healthcare & Medical Device

In the healthcare sector, there has been an expansion in the use of 3D printing in a number of different applications. For some time now, 3D printing has been leveraged in the production of customized prosthetics and dental and orthopedic implants, as well as hearing aids. People no longer have to wait months for a needed prosthetic, which is particularly important for children, who in the past could outgrow their prosthetic before receiving it.





3D-printed implants are also fairly commonplace; for example, more than 100,000 hip replacements produced by GE Additive have used this technology. In addition to dental implants, molds for clear aligners are among the objects most frequently produced via 3D printing around the globe. In addition, most hearing aids produced in the last 20 years have also been made using additive manufacturing. 3D-printed models of patient organs are also used by surgeons for practice before performing complicated operations.

### Pharmaceutical

Pharmaceutical drug development and manufacturing may also be impacted by the introduction of 3D-printing technology. The ability to produce personalized, 3D-printed drug products on demand would help speed drug development by providing rapid access to different dosage forms for preclinical, first-in-human clinical studies and later-stage clinical trials. It would also make highly personalized medications accessible to patients through dramatically streamlined customized manufacturing, an approach that is simply not possible today.

In early clinical trials, for instance, dose-escalating studies could be expedited with the ability to rapidly 3D print different formulated products with different and assured API loadings and release kinetics. Similarly, once a drug was approved, 3D printing would make it possible to individualize drug products for different patients and adjust those products with varied quantities of drug substance to meet changing patient needs.

Indeed, 3D printing pharmaceutical tablets is attractive because it can leverage well-established hot-melt extrusion (HME) technology, which is increasingly used to produce drug products with improved solubility and bioavailability properties. Instead of compressing the generated granules into tablets or pellets to fill capsules, they could be formed into filaments containing specific API loadings and/or excipient contents for use in 3D printers.

MilliporeSigma is exploring an even more efficient approach that leverages melt drop deposition as the 3D printing technology. This method allows the direct use of the HME granules without the need to form filaments.

So far, one 3D printed drug — Spritam (levetiracetam, Aprexia Pharmaceuticals) — has received FDA approval, and a second — T19 (Triastek) — received FDA approval for its Investigational New Drug (IND) application in February 2021.

Ongoing work at major and emerging pharma companies is largely focused on combining additive manufacturing with other advanced digital solutions, from artificial intelligence to formulation automation, to further optimize formulations for 3D printing. Initial applications will likely be for small-volume therapies targeting small patient populations, but given the scalable nature of additive manufacturing, there is potential for its use in the decentralized production of tableted drug products, and perhaps even on-demand production in clinics or at home.

It is also worth noting that [3D printing technology of biological materials](#) (e.g., cells, media, nutrients, cytokines) is rapidly advancing and has applications ranging from tissue engineering to drug discovery and disease diagnosis. More specifically, bioprinting is being explored for the production of replacement organs to avoid the need for organ donation and for the printing of organ models that can be used instead of animals for drug screening during development and disease diagnosis for patients. This technology is also being investigated as an alternative means for meat production.

## Method and Material Impacts on Cleaning Requirements

Each material used in 3D printing comes with its own set of strengths and limitations, with many requiring post-processing steps to finalize the print. In addition, while many parts of 3D printers are similar and have similar cleaning requirements, each 3D-printing method has unique critical cleaning requirements. Alconox Inc. offers solutions that can help your prints become as clean as possible when finished and keep 3D printers clean and performing at a high level.

### Materials

Extruded plastics —like ABS (acrylonitrile butadiene styrene), PLA (polylactic acid), and nylon — are among the most common materials used in household printers as well as industrial settings. A scrub or soak in a 1–2% solution of a mild, free-rinsing, alkaline detergent to remove residual adhesive, small deposits of material, stringing, and any dirt or dust, will provide a critically clean print after printing or removal of the supports. A wash after sanding to remove any particulates is also recommended. As some plastics can be heat sensitive when combined with aqueous detergents, it is further recommended that wash temperatures be lukewarm until compatibility is confirmed.

After each of the two curing processes with photosensitive resins, excess liquid resin must be removed from the print. Cleaning can be achieved with either an aqueous solution or in a solvent bath, often isopropyl alcohol (IPA).






Using solvents like IPA comes with a number of challenges, however, including flammability and waste concerns, cost, and effectiveness due to only a single solvent-based cleaning mechanism. A 1–2% solution of a mild alkaline aqueous detergent, on the other hand, is able to clean off most liquid resins without harming the printed piece, and without any hazardous waste generation or flammability concerns.

Metals and metal alloy-based 3D printed objects are produced in a bed of powder and often have a layer of powder adhering to their surfaces when pulled from the bed. A detergent bath featuring some agitation can be a great help in removing this layer. Avoiding corrosion is important with metals, so picking a compatible detergent, and quick removal of rinse water for reactive metals is essential.

Dissolvable supports are used to hold up complex prints without leaving unsightly marks and require water or other solvents to remove. Submerging the entire printed part/support assembly in a detergent bath can often help with removal of these supports and any residue the dissolved material may leave on the print.



## Printers

Printer beds in 3D printers are typically composed of polyetherimide (PEI) sheets, glass, or an adhesive tape (on non-heated beds). Cleaning them on a regular basis will minimize the chances of any imperfections.

Cleaning of printer beds is similar for all types of printers. Glass beds and beds with a PEI sheet can be cleaned using IPA or other solvents, such as acetone. Acetone, like IPA, will have the same challenges, including flammability and waste concerns, cost, and effectiveness due to only a single solvent-based cleaning mechanism. A mild alkaline, aqueous detergent would again be the safer, more effective option.

Any remaining adhesive should then be washed and gently scraped off. The bed can then be washed with a detergent solution to remove any remaining residue. Older PEI sheets that have seen a lot of use may also benefit from sanding with a fine-grit (1,500+) sandpaper. Beds with adhesive tape are the most difficult to clean, as the tape is subject to breakage. Gentle scraping with a spatula is most effective. These processes may also lend to the generation of particulate residue that mild alkaline aqueous detergents can effectively address.

With vat polymerization printers, the main components that require cleaning are the light-emitting/focusing system and the resin vat. Cleaning the laser emitter and/or mirrors can be performed with a clean cloth and some mild alkaline detergent solution followed by a light rinse to ensure that no dirt or residue is present. The resin vat should be thoroughly cleaned when changing types of resin so that you don't get mixed resins and cross contamination of materials.

For powder-bed printers, it is important to keep the print head or laser system and the powder bed itself clean. Cleaning the print head is important for proper deposition of the binding agent or dye. It is recommended to inspect the print head for obstructions and gently clean with a stiff brush and mild alkaline detergent. Cleaning out the powder bed is generally only a concern when changing materials, as it is undesirable to have leftover powder mixed in with the new, leading to cross contamination.



## Printer Nozzles

Fused deposition modeling printers that operate via extrusion also require cleaning of the extruder nozzle to prevent blockages. The exterior of the nozzle can be cleaned while still hot using a damp cloth with detergent if necessary. If that is not sufficient, a wire brush, small blade, or needle can be used.

To prevent nozzle malfunction and blockage, you first need to make sure that the tip of the nozzle and the area around the heating block is clean

Signs of a blocked nozzle include:

- The filament is not extruding uniformly.
- The nozzle extrudes very thin filament.
- Nothing comes out from the nozzle.

Before you begin, you'll need a 2% solution of a mild alkaline detergent, a torch, and a very thin wire. To remove and clean a blocked extruder nozzle:

1. Soak the removed nozzle in the warm detergent for about 15 minutes to clean out exterior dirt. Rinse thoroughly. Use a soft cloth to clean the nozzle.
2. Place the nozzle on a stone and burn it using the torch for about one minute. Make sure it is extremely hot. You should see slight changes in the color.
3. Use a very thin wire to clear the hole in the nozzle. If the wire cannot go through, repeat step 2 until it can go through. Don't force through the hole with the wire.

4. If the jam was not caused by ABS or another filament that will dissolve or be cleaned away, you could also use a heat gun or blowtorch to heat it and melt the clog away. Once the nozzle is very hot, you can use a thin, soft wire to check if you can clear a hole through the nozzle.

### Support Material


For cleaning soluble support material, you can use a selected choice of alkaline detergents in an ultrasonic machine. For ultrasonic cleaning of all the resin from around and inside your 3D prints, we recommend a 2% [ALCONOX® Powdered Precision Cleaner](#) or [DETNOX® Ultimate Precision Cleaner](#); warm temperatures (>130°F/55°C); and an ultrasonic bath. Our powerful emulsifying detergent solutions, along with sonication, would be most effective in removing the residue, especially particulates.

For a phosphate-free alternative, [LIQUINOX® Critical Cleaning Liquid Detergent](#) is an optimal choice.

All three of the above detergents are compatible for cleaning a wide range of polymeric 3D-printing substrates and helping with complete removal of soluble support materials. (It is highly recommended that you verify the cleaning approach in your own setup, conditions, application, and instruments.)

Some sensitive plastics are susceptible to [stress cracking](#), in which case a surfactant-free detergent would be needed. We would highly recommend [DETOJET® Low Foaming Liquid Detergent](#) in such a case. Typically, sensitive or stressed plastics are not used in 3D printing applications.

In terms of safety, it is highly recommended to avoid using flammable liquids like alcohol and acetone in your ultrasonic cleaner. There is a risk of your ultrasonic cleaner causing a small spark, and that would be enough to cause a fire. If you have an ultrasonic transducer that fails, the energy from it can transfer into the cleaning fluid, which if flammable, this can also result in a fire.



### Cleaning Chemistry

Since most 3D printing-related soils are light oils and particulates, mild alkaline cleaners are very effective. Alkaline cleaners remove organics, including oils, resins, extracts, and an array of other soils. Aqueous, mild alkaline detergents containing micelle-forming surfactants for emulsifying provide a safe, efficient, effective choice.

Alconox Inc. is armed with 75 years of experience recommending cleaners for use in [laboratory](#), [cosmetics](#), [medical devices](#), [food](#), [biotech](#), and [pharmaceuticals](#), where we have gained experience cleaning an array of residues and resins. This same experience allows us to bring great insight into removing residues found in the additive manufacturing process.

When large amounts of residue are attached to manufacturing surfaces, increased detergent concentration is generally required. Typically, a 1–2% (10–20 mL/L or 1.3–2.5 oz/gal) concentration suffices for routine cleaning; however, 2–5% (20–50 mL/L or 2.5–6 oz/gal) is often required for stickier or more dense residues.

Heat is an important variable in the cleaning of photosensitive resin prints. Within reason, higher temperatures increase the mobility of resins via softening, and accelerates the cleaning process residue removal. Using the highest temperature compatible with your prints to be cleaned ensures efficient emulsification of any leftover resin.

### Cleaning Methods

Cleaning methods can be generally divided into manual or clean out of place (COP) cleaning, and automated cleaning.

Manual and COP cleaning includes methods that include direct, applied mechanical action. This is where higher-foaming detergents are not only permissible but useful, as they are best for emulsifying. Emulsifying is the preferred approach to remove organic soils, like oils, extracts, and resins. Manual/COP cleaning methods are soaking, scrubbing, and sonication (ultrasonic tanks).

Automated methods include the use of washers, high-pressure sprayers, and clean-in-place (CIP) systems. Low-foaming detergents use surfactants, which foam less and thereby eliminate situations where over-foam can occur. The less efficient surfactants are compensated for by higher alkalinity, chelation, wetting agents, and other methods to ensure that residues are quickly, safely, and reliably removed.

### Manual Cleaning / COP Detergents

We would recommend the following detergents for soaking, scrubbing, or sonication of your parts from particulates, resins, and oily residues:

#### [ALCONOX® Powdered Precision Cleaner](#)

is a powerful emulsifier of organic and oily residues, used for many decades in manufacturing cleaning. It is biodegradable, drain-safe, and completely free-rinsing. An excellent choice for general cleaning of manual and COP 3D printing applications.





### **DETONOX® Ultimate Precision Cleaner**

is our most potent manual detergent. A liquid concentrate, ideal and safe for hand and ultrasonic use on difficult 3D printing–related residues. It is a non-caustic detergent for exceptional removal of difficult, adherent resin and particulates.

### **LIQUINOX® Critical Cleaning Liquid Detergent**

is a phosphate-free liquid alternative to Alconox detergents for manual and COP cleaning.

### **CITRANOX® Liquid Acid Cleaner and Detergent**

is a phosphate-free acidic cleaner and metal brightener for the removal of oxidation, scale, salts, and inorganic residues in manual and ultrasonic cleaning. Effective for removal of 3D printing residues in soft and reactive metals.

For each of the approaches above, we would recommend warm to high-heat temperatures (140°F / 60°C+) and a 1–2% starting concentration. Initial rinse should also be of a similar temperature to avoid thermal shock of the micelles that form.

### **CIP System and Washer Detergents**

For cleaning via washers and other automated high-pressure units (e.g., CIP) necessitating low-foaming detergents, we would recommend:



### **ALCOJET® Low Foaming Powdered Detergent**

is ideal for removing 3D printing residues in washers with a cup-in-door designed to hold powdered detergent, often with a door that stays closed during the prewash cycle and then pops open during the wash cycle.

### **SOLUJET® Low-Foaming Phosphate-Free Liquid**

is a low foaming liquid concentrate detergent for removal of 3D printing residues in spray CIP system and washers with a liquid dispensing system.

### **CITRAJET® Low-Foam Liquid Acid Cleaner/Rinse**

is a low-foaming, phosphate-free acidic cleaner and metal brightener for cleaning of 3D printing residues from sensitive and reactive metals in low foaming applications including CIP and washer systems.

### **DETOJET® Low Foaming Liquid Detergent**

is ideal for cleaning sensitive plastics and stressed polycarbonates with 3D printing residues, as it contains no surfactants and will help avoid stress cracking.

A concentration of 1–2% detergent (10–30 mL/L or 1.25–2 oz/gal) is a good place to start, and the highest temperature that can be reasonably achieved. Recall that heat expedites cleaning and higher concentrations increase capacity or the amount that can be removed by a detergent. Therefore, qualitatively speaking, if there is a large amount of residue, higher concentrations may be needed.

**For more information, see [Critical Cleaning for 3D Printing in Precision Manufacturing](#).**





## Alconox Inc. Product Innovation

**Alconox Inc. is honored to have reached our 75th anniversary milestone. Lets take a look back at some important dates in the company's history.**

**Pre 1941** While working at a hosiery company, chemist Louis Zisman learned about the emerging chemistry of synthetic surfactants. He discovered that a wetting agent he developed excelled at cleaning glassware without leaving a residue.

**1940s** Louis Zisman combined his talents with that of William Lebowitz, owner of Standard Scientific Supply Corporation. Standard Scientific was a distributor of products focused on the laboratory market. Together, they began selling and distributing detergent to hospitals and laboratories.

The detergent was first used in commerce on September 1st, 1943 and was coined, ALCONOX® Powdered Precision Cleaner.

Alconox Inc. was formed May 28th, 1946, to hold exclusive rights to ALCONOX detergent and its distribution.

**1950s** Alconox Inc. opens an office on 853 Broadway, New York, NY. Over time, more and more dealers signed up to distribute Alconox detergent and the product became widely distributed by just about every lab or hospital dealer.

Alconox Inc. launches a decade of innovation with some hits and some misses. They created the powdered detergent, ORANOX® Rust Inhibitor, that targeted dentists and oral surgeons. It was ultimately not successful and was discontinued in 1954. Recently, one box was recovered from an antique store in Albuquerque, New Mexico.

ALCOTABS® Tablet Pipet Detergent, a tablet solution to better clean siphon pipet washers, was developed and hit the market July 12th, 1957. The new, innovative, flotation cleaning, slow dissolving tablet allowed for cleaning over multiple cycles in the pipet washer resulting in better pipet cleaning.

Closing out the decade came another Alconox Inc. powerhouse, in fact our flagship detergent for lab washers. ALCOJET® Low Foaming Powdered Detergent was created and put in distribution August 12th, 1959. It was invented to provide a low-foaming detergent effective in automatic dishwasher glassware washing applications.

**1960s** On Oct. 28th, 1964 LIQUINOX® Critical Cleaning Liquid Detergent, was developed in response to customer requests for a phosphate-free detergent that would not interfere with biological oxygen demand (BOD) bottle washing where unrinsed bottles could result in contamination with phosphate essential nutrients causing false high readings in BOD testing.

**1970s** In order to improve removal of protein residues of blood, body fluids, and tissue for cleaning medical instruments, TERGAZYME® Enzyme-Active Powdered Detergent was introduced on November 21st, 1971.

In order to improve removal of protein residues of blood, body fluids, and tissue for cleaning medical instruments, TERGAZYME® Enzyme-Active Powdered Detergent was introduced on November 21st, 1971.

**1980s** In response to requests for a nonionic and ion-free detergent that does not leave conductive residues on electronic components, DETERGENT 8® Low-Foaming Ion-Free Detergent was introduced on July 14th, 1983.

For improved inorganic salts, oxides and acid labile residue cleaning, CITRANOX® Liquid Acid Cleaner and Detergent was introduced on December 8th, 1987.

**1990s** To help replace hazardous ozone depleting solvents that were being restricted by the Montreal Protocol, LUMINOX® Low-Foaming Neutral pH Liquid Detergent was first sold on April 21st, 1997.

**2000s** Per the need for a low foaming acid detergent and neutralizing rinse for automatic lab washers, CITRAJET® Low-Foam Liquid Acid Cleaner/Rinse was released on March 13th, 2001.

To address requests for an automatic lab washer phosphate free detergent for use in cleaning biological oxygen bottles (BOD), SOLUJET® Low-Foaming Phosphate-Free Liquid, went on sale July 2004.

Concurrently, TERGAJET® Low-Foaming Phosphate-Free Detergent was also released at the same time in order to have a powdered product for use in automated washer BOD bottle cleaning.

TERGAZYME makes a cameo on *Bones* Season 10 Episode 4.

Working with customers with difficult cosmetic polymer and sticky botanical extract residues we developed DETONOX® Ultimate Precision Cleaner, introduced in April 2015.

To match the manual cleaning power of DETONOX for cosmetic polymers and sticky botanical residues in an automated washer cleaning, KEYLAJET® Low-foaming Chelating High Alkaline Liquid was released on August 18th, 2016.

Alconox Inc. celebrates 75 years of solving critical cleaning challenges.



Critical Cleaning Experts

30 Glenn Street, Suite 309  
White Plains, NY 10603 • USA

## Get Validation Support or Help With Your Critical Cleaning Challenge

Alconox Inc. has more than 75 years' experience developing aqueous cleaning solutions for pharmaceutical manufacturing. Let us help solve your next critical cleaning challenge.

Please contact Alconox Inc. for expert validation support or verification laboratory services:

[cleaning@alconox.com](mailto:cleaning@alconox.com)



## Learn More About Critical Cleaning

Request a FREE copy of:

**The Aqueous Cleaning Handbook**  
or  
**Critical Cleaning Guide**



## Try a Free Sample of Alconox Inc. Detergents

Use our sample request form at [alconox.com](http://alconox.com). Or call:

**++914-948-4040**

**For questions or comments about this white paper, please contact Alconox Inc. Technical Support at 914.948.4040 or [cleaning@alconox.com](mailto:cleaning@alconox.com)**