



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/12/2019

## ***CERTIFICATE OF ANALYSIS*** **ALCONOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: A1D9 A1D9A, A1D9B, A1D9C, A1D9D, A1D9E

Sub lot numbers of Master lot No: A2H8

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 4/11/2019

Expires end of: 04/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.5	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	32.1	29.0 - 33.0

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*