



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/4/2016

## **CERTIFICATE OF ANALYSIS**

### **ALCONOX**

#### **Manufactured by: Alconox, Inc.**

Sub lot No: A1E6 A1E6A, A1E6B, A1E6C, A1E6D, A1E6E

Sub lot numbers of Master lot No: A2M5

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 5/2/2016

Expires end of: 05/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

#### **ANALYSIS**

#### **RESULT**

#### **ACCEPTED RANGE**

Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.2	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	30.7	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*