



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/1/2018

CERTIFICATE OF ANALYSIS
ALCONOX

Manufactured by: Alconox, Inc.

Sub lot No: A5J8 A5J8A, A5J8B, A5J8C, A5J8D, A5J8E

Sub lot numbers of Master lot No: A3L7

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 9/27/2018

Expires end of: 09/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.4	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	31.4	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.