



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/21/2017

CERTIFICATE OF ANALYSIS ***ALCOTABS***

Manufactured by: Alconox, Inc.

Sub lot No: B1D7

Sub lot numbers of Master lot No: B1D6

Manuf. Cat. no.: 1500, 1500-1

Date Of Manufacture: 4/1/2017

Expires end of: 04/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White and cream colored tablet
pH: (1% solution in DI water)	6.8	5.9 to 8.8
Surface Tension:(1% solution Dyne/cm)	30.7	29.0 - 33.0
Dissolve Time	16.4	Effervescent in water for a minimum of 7 minutes

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.