



Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave, Jersey City 07304 USA

04/30/2018

CERTIFICATE OF ANALYSIS
20% CITRIC ACID

Manufactured by: Alconox, Inc.

Lot No: C041618

Manuf. Cat. no.: C201, C201-1

Date of Manufacture: 04/27/2018

Shelf Life Expiration Date: 04/2020

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Analysis</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear liquid
Refractive Index (degrees brix)	18.5	16.5-21.0
Specific Gravity:	1.094	1.070-1.100

Malcolm McLaughlin

Alconox, Inc