



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/23/2018

## CERTIFICATE OF ANALYSIS CITRANOX

**Manufactured by: Alconox, Inc.**

Sub lot No: C1B8 C1B8F, C1B8J, C1B8K, C1B8L, C1B8M

Sub lot numbers of Master lot No: C2L6

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 2/20/2018

Expires end of: 02/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u>                       | <u>RESULT</u> | <u>ACCEPTED RANGE</u>       |
|---------------------------------------|---------------|-----------------------------|
| Appearance:                           | OK            | Clear, olive colored liquid |
| pH: (1% solution in DI water)         | 2.8           | 2.0 - 3.0                   |
| Surface Tension:(1% solution Dyne/cm) | 31.6          | 29.0 - 33.0                 |
| Specific Gravity:                     | 1.116         | 1.09 - 1.14                 |

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*