



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/7/2017

## **CERTIFICATE OF ANALYSIS** **CITRANOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: C1D7 C1D7F, C1D7J, C1D7K, C1D7L, C1D7M

Sub lot numbers of Master lot No: C1A5

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 4/6/2017

Expires end of: 04/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, olive colored liquid
pH: (1% solution in DI water)	2.9	2.0 - 3.0
Surface Tension:(1% solution Dyne/cm)	31.3	29.0 - 33.0
Specific Gravity:	1.114	1.09 - 1.14

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*