



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/8/2018

CERTIFICATE OF ANALYSIS **CITRANOX**

Manufactured by: Alconox, Inc.

Sub lot No: C1F8 C1F8F, C1F8J, C1F8K, C1F8L, C1F8M

Sub lot numbers of Master lot No: C2K7

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 6/7/2018

Expires end of: 06/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---------------------------------------|-------|-----------------------------|
| Appearance: | OK | Clear, olive colored liquid |
| pH: (1% solution in DI water) | 2.6 | 2.0 - 3.0 |
| Surface Tension:(1% solution Dyne/cm) | 30.3 | 29.0 - 33.0 |
| Specific Gravity: | 1.112 | 1.09 - 1.14 |

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.