



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/26/2019

**CERTIFICATE OF ANALYSIS**  
**CITRANOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: C1F9 C1F9F, C1F9J, C1F9K, C1F9L, C1F9M

Sub lot numbers of Master lot No: C1F8

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 6/25/2019

Expires end of: 06/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, olive colored liquid
pH: (1% solution in DI water)	2.7	2.0 - 3.0
Surface Tension:(1% solution Dyne/cm)	31.8	29.0 - 33.0
Specific Gravity:	1.11	1.09 - 1.14

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*