



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/1/2018

**CERTIFICATE OF ANALYSIS**  
**CITRANOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: C1J8 C1J8F, C1J8J, C1J8K, C1J8L, C1J8M

Sub lot numbers of Master lot No: C2K7

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 9/28/2018

Expires end of: 09/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u>                          | <u>RESULT</u> | <u>ACCEPTED RANGE</u>       |
|--|---------------|-----------------------------|
| Appearance:                              | OK            | Clear, olive colored liquid |
| pH: (1% solution in DI water)            | 2.9           | 2.0 - 3.0                   |
| Surface Tension:(1% solution<br>Dyne/cm) | 31.8          | 29.0 - 33.0                 |
| Specific Gravity:                        | 1.115         | 1.09 - 1.14                 |

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*