



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/4/2019

CERTIFICATE OF ANALYSIS
CITRANOX

Manufactured by: Alconox, Inc.

Sub lot No: C2A9 C2A9F, C2A9J, C2A9K, C2A9L, C2A9M

Sub lot numbers of Master lot No: C2K7

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1830, 1855

Date Of Manufacture: 1/30/2019

Expires end of: 01/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|--|------|-----------------------------|
| Appearance: | OK | Clear, olive colored liquid |
| pH: (1% solution in DI water) | 2.9 | 2.0 - 3.0 |
| Surface Tension:(1% solution Dyne/cm) | 31.4 | 29.0 - 33.0 |
| Specific Gravity: | 1.12 | 1.09 - 1.14 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.