



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date Of Manufacture: 05/20/2016
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C9E6
BCI Lot number: 622831
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 04/30/2018 To meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.8	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	31.4	29.0 – 33.0
Specific Gravity	1.112	1.090 – 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	3	Not Applicable
Number of Bar Code Labels Discarded	14	Not Applicable


 Malcolm McLaughlin, QA

REF 81912

 2018-04-30 **LOT** 622831



+H628819121A



+\$2043018622831A/



(01)15099590223571
 (17)180430
 (10)622831

Alconox, Inc.
The Leader in Critical Cleaning
 30 Glenn St. White Plains NY 10603 USA
 Tel 914.948.4040 • Fax 914-948-4088
 www.alconox.com • cleaning@alconox.com