



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date Of Manufacture: 07/17/2017
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C9E7
BCI Lot number: 123162
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 07/31/19 to meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.9	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	31.8	29.0 - 33.0
Specific Gravity	1.114	1.090 - 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	4	Not Applicable
Number of Bar Code Labels Destroyed	24	Not Applicable

Malcolm McLaughlin, QA

REF 81912

2019-07-31


LOT 123162



H628819121A



S52073119123162A



(01)15099598223571
(17)190731
(10)123162

Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com