



**CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)**

Date Of Manufacture: 04/10/2018
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C1D8
BCI Lot number: C9F8
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 03/31/20 to meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>-</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.8	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		31.8	29.0 - 33.0
Specific Gravity		1.108	1.090 - 1.140
Number of Bar Code Labels Used for Labels		720	Not Applicable
Number of Bar Code Labels Used for Records		3	Not Applicable
Number of Bar Code Labels Destroyed		519	Not Applicable


Malcolm McLaughlin, QA

REF 81912

 2020-03-31 **LOT** C9F8



H628019221A



S52033120C9F8AC



(01)15099590223571
(17)200331
(10)C9F8

Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com