



**CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)**

Date of Manufacture: 05/20/2018
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C1E8
BCI Lot number: C9G8
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 05/31/20 to meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.7	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	31.9	29.0 - 33.0
Specific Gravity	1.118	1.090 - 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	3	Not Applicable
Number of Bar Code Labels Destroyed	431	Not Applicable

Malcolm McLaughlin, QA

REF 81912

2020-05-31

LOT C9G8



+H628819121A



+52053120C9G8AF



(01)15099590223571
(17)200531
(10)C9G8

Alconox, Inc.

The Leader in Critical Cleaning

30 Glenn St. White Plains NY 10603 USA

Tel 914.948.4040 • Fax 914-948-4088

www.alconox.com • cleaning@alconox.com