



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/16/2016

***CERTIFICATE OF ANALYSIS***  
***DETOJET***

**Manufactured by: Alconox, Inc.**

Sub lot No: D1F6 D1F6F, D1F6G, D1F6H, D1F6J, D1F6K, D1F6L, D1F6M

Sub lot numbers of Master lot No: D1A5

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 6/10/2016

Expires end of: 06/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.3	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.338	1.305 - 1.358

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*