



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 12/02/2011

## **CERTIFICATE OF ANALYSIS** **DETERGENT 8**

**Manufactured by: Alconox, Inc.**

Sub lot No: E1F9 E1F9G, E1F9J, E1F9K, E1F9L, E1F9M

Sub lot numbers of Master lot No: E1C8

Manuf. Cat. no.: 1701, 1701-1, 1705, 1715, 1755

Date Of Manufacture: 6/11/2009

Shelf Life Expiration Date (1095 days): 6/10/2012

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	11.0	10.5 - 11.8
Surface Tension:(1% solution Dyne/cm)	32.1	30.0 - 34.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	0.992	0.970 - 1.010

**Malcolm McLaughlin**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*