



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 12/19/2017

CERTIFICATE OF ANALYSIS **TERGAZYME**

Manufactured by: Alconox, Inc.

Sub lot No: G1M7 G1M7A, G1M7B, G1M7C, G1M7D

Sub lot numbers of Master lot No: G1L6

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 12/18/2017

Expires end of: 12/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White powder with cream colored flakes and small brown granules
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.4	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	30.7	29.0 - 33.0

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.