



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 1/18/2018

## **CERTIFICATE OF ANALYSIS** **ALCOJET**

**Manufactured by: Alconox, Inc.**

Sub lot No: J1A8 J1A8A, J1A8B, J1A8C, J1A8D

Sub lot numbers of Master lot No: J1K6

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 1/16/2018

Expires end of: 01/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.1	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	32.1	31.0 - 36.0
Foam (1% in DI water)	NIL	NIL

**Malcolm McLaughlin, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*