



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/26/2019

CERTIFICATE OF ANALYSIS
ALCOJET

Manufactured by: Alconox, Inc.

Sub lot No: J1G9 J1G9A, J1G9B, J1G9C, J1G9D

Sub lot numbers of Master lot No: J1L8

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 7/23/2019

Expires end of: 07/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u> | <u>RESULT</u> | <u>ACCEPTED RANGE</u> |
|---|----------------------|------------------------------|
| Appearance: | OK | White powder |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 11.5 | 11.0 - 12.5 |
| Surface Tension:(1% solution Dyne/cm) | 32.5 | 31.0 - 36.0 |
| Foam (1% in DI water) | NIL | NIL |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.