



Alconox, Inc.

The Leader in Critical-Cleaning  
30 Glenn St., Suite 309, White Plains, NY 10603 USA  
Tel.914.948.4040 • Fax.914.948.4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/17/2016

## **CERTIFICATE OF ANALYSIS**

### **ALCOJET**

**Manufactured by: Alconox, Inc.**

Sub lot No: J1K6 J1K6A, J1K6B, J1K6C, J1K6D

Sub lot numbers of Master lot No: J1L5

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 10/13/2016

Expires end of: 10/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.5	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	32.3	31.0 - 36.0
Foam (1% in DI water)	NIL	NIL

**Malcolm McLaughlin, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*