



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/23/2018

CERTIFICATE OF ANALYSIS **ALCOJET**

Manufactured by: Alconox, Inc.

Sub lot No: J1K8 J1K8A, J1K8B, J1K8C, J1K8D

Sub lot numbers of Master lot No: J1K7

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 10/18/2018

Expires end of: 10/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.2	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	33.3	31.0 - 36.0
Foam (1% in DI water)	NIL	NIL

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.