



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 11/8/2019

***CERTIFICATE OF ANALYSIS***  
***ALCOJET***

**Manufactured by: Alconox, Inc.**

Sub lot No: J1L9 J1L9A, J1L9B, J1L9C, J1L9D

Sub lot numbers of Master lot No: J1L8

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 11/4/2019

Expires end of: 11/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.6	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	32.9	31.0 - 36.0
Foam (1% in DI water)	NIL	NIL

**Laura Falla, QA**

**Alconox, Inc.**

***Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.***