



Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave, Jersey City 07304 USA

07/29/2016

CERTIFICATE OF ANALYSIS
Keylajet

Manufactured by: Alconox, Inc.

Lot No: K1G6

Manuf. Cat. no.: 2405, 2455

Date of Manufacture: 07/25/2016

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.

Expires End of: 07/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>Analysis</u>	<u>Result</u>	<u>Acceptance Range*</u>
Appearance:	OK	Clear to yellow liquid
pH: (1% solution in DI water)	12.76	12.65 +/- 0.8
Specific Gravity	1.368	1.37 +/- 0.05

*Acceptance ranges are based on experience and subject to change from statistical process control

Malcolm McLaughlin

Alconox, Inc, QA