



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/25/2018

CERTIFICATE OF ANALYSIS

SLQ Stryker - Liquinox

Manufactured by: Alconox, Inc.

Sub lot No: L022618

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: 9701, 9705, 9715, 9755

Date Of Manufacture: 2/26/2018

Expires end of: 02/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear golden liquid
pH: (as is)	8.3	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.7	29.0 - 33.0
Specific Gravity:	1.080	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.