



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/11/2016

## ***CERTIFICATE OF ANALYSIS*** ***LIQUINOX***

**Manufactured by: Alconox, Inc.**

Sub lot No: L3C6 L3C6F, L3C6H, L3C6J, L3C6K, L3C6L, L3C6M

Sub lot numbers of Master lot No: L2A5

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 3/28/2016

Expires end of: 03/18

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.5	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.1	29.0 - 33.0
Specific Gravity:	1.078	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*