



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/26/2019

CERTIFICATE OF ANALYSIS ***LIQUINOX***

Manufactured by: Alconox, Inc.

Sub lot No: L3G9 L3G9F, L3G9H, L3G9J, L3G9K, L3G9L, L3G9M

Sub lot numbers of Master lot No: L1J8

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 7/25/2019

Expires end of: 07/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.7	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	30.4	27.0 - 31.0
Specific Gravity:	1.07	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.