



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 8/27/2019

## ***CERTIFICATE OF ANALYSIS*** ***LIQUINOX***

**Manufactured by: Alconox, Inc.**

Sub lot No: L3H9 L3H9F, L3H9H, L3H9J, L3H9K, L3H9L, L3H9M

Sub lot numbers of Master lot No: L1J8

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 8/26/2019

Expires end of: 08/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.3	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	30.4	27.0 - 31.0
Specific Gravity:	1.07	1.06 - 1.09

**Laura Falla, QA**

**Alconox, Inc.**

***Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.***