



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/1/2018

**CERTIFICATE OF ANALYSIS**  
**LIQUINOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: L4J8 L4J8F, L4J8H, L4J8J, L4J8K, L4J8L, L4J8M

Sub lot numbers of Master lot No: L1H7

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 9/28/2018

Expires end of: 09/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.6	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	30.2	27.0 - 31.0
Specific Gravity:	1.075	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*