



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 8/22/2017

**CERTIFICATE OF ANALYSIS**  
**SLQ - Stryker Liquinox**

**Manufactured by: Alconox, Inc.**

Lot No: L8H7

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: 9715, 9755

Date Of Manufacture: 8/18/2017

Expires end of: 08/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.3	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.6	29.0 - 33.0
Specific Gravity:	1.076	1.06 - 1.09

Malcolm McLaughlin, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*