



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 9/3/2019

**CERTIFICATE OF ANALYSIS**  
**LUMINOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: M1H9 M1H9G, M1H9J, M1H9L, M1H9M

Sub lot numbers of Master lot No: M1L8

Manuf. Cat. no.: 1901, 1901-1, 1915, 1955

Date Of Manufacture: 8/16/2019

Expires end of: 08/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear liquid
pH: (as is)	6.9	6.0 - 8.0
Surface Tension:(1% solution Dyne/cm)	35.7	32.0 - 38.0
Specific Gravity:	1.02	1.00 - 1.06
Foam (1% in DI water)	NIL	NIL

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*