



Alconox, Inc.  
The Leader in Critical Cleaning  
30 Glenn St. White Plains NY 10603 USA  
Tel 914.948.4040 • Fax 914-948-4088  
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave, Jersey City 07304 USA

2/1/2016

**CERTIFICATE OF ANALYSIS**  
**LQ9**

Manufactured by: Alconox, Inc.

Lot No: Q2A6

Manuf. Cat. no.: LQ941, LQ905, LQ915, LQ930, and LQ955

Date of Manufacture: 1/22/2016

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*

Shelf Life Expiration Date: 1/21/2018

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Analysis</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear golden liquid
pH: (as is)	8.4	8.0-9.5
Surface Tension: (1% solution dyne/cm)	31.4	29.0 – 33.0
Specific Gravity:	1.078	1.060 – 1.090

Malcolm McLaughlin

Alconox, Inc