



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/8/2019

## CERTIFICATE OF ANALYSIS SPCC5

**Manufactured by: Alconox, Inc.**

Sub lot No: S020419

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: SPCC5-5, SPCC5-55

Date Of Manufacture: 2/4/2019

Expires end of: 02/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

|  |      |   |
|--|------|---|
| Appearance:                              | OK   | Clear pale yellow liquid –nearly odorless |
| pH: (2.5% solution in DI water)          | 10.2 | 9.5 - 11.5                                |
| Surface Tension:(1% solution<br>Dyne/cm) | 30.2 | 29.0 - 33.0                               |
| Specific Gravity:                        | 1.18 | 1.155 -1.225                              |

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*