



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/25/2017

CERTIFICATE OF ANALYSIS

SPCC5

Manufactured by: Alconox, Inc.

Sub lot No: S042517

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: SPCC5-5, SPCC5-55

Date Of Manufacture: 4/25/2017

Expires end of: 04/19

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear pale yellow liquid –nearly odorless
pH: (2.5% solution in DI water)	10.2	9.5 - 11.5
Surface Tension:(1% solution Dyne/cm)	29.3	29.0 - 33.0
Specific Gravity:	1.174	1.155 -1.225

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.