



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 - Fax.914.948.4088

www.alconox.com - cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 11/26/2019

***CERTIFICATE OF ANALYSIS***  
***SPCC5***

**Manufactured by: Alconox, Inc.**

Sub lot No: S112519

Sub lot numbers of Master lot No: None

Manuf. Cat. no.: SPCC5-5, SPCC5-55

Date Of Manufacture: 11/25/2019

Expires end of: 11/21

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear pale yellow liquid –nearly odorless
pH: (2.5% solution in DI water)	10.1	9.5 - 11.5
Surface Tension:(1% solution Dyne/cm)	30.8	29.0 - 33.0
Specific Gravity:	1.18	1.155 -1.225

**Laura Falla, QA**

**Alconox, Inc.**

***Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.***