



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/12/2018

CERTIFICATE OF ANALYSIS **SOLUJET**

Manufactured by: Alconox, Inc.

Sub lot No: S1G8 S1G8F, S1G8G, S1G8J, S1G8K, S1G8L, S1G8M

Sub lot numbers of Master lot No: S1J7

Manuf. Cat. no.: 2101, 2101-1, 2105, 2115, 2130, 2155

Date Of Manufacture: 7/11/2018

Expires end of: 07/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	pale straw colored liquid
pH: (1% solution in DI water)	12.2	11.5 - 13.0
Surface Tension:(1% solution Dyne/cm)	32.3	31.5 - 36.5
Specific Gravity:	1.280	1.26 - 1.34

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.