



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 4/27/2018

## **CERTIFICATE OF ANALYSIS** **SOLUJET**

**Manufactured by: Alconox, Inc.**

Sub lot No: S2D8 S2D8F, S2D8G, S2D8J, S2D8K, S2D8L, S2D8M

Sub lot numbers of Master lot No: S1J7

Manuf. Cat. no.: 2101, 2101-1, 2105, 2115, 2130, 2155

Date Of Manufacture: 4/26/2018

Expires end of: 04/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	pale straw colored liquid
pH: (1% solution in DI water)	12.1	11.5 - 13.0
Surface Tension:(1% solution Dyne/cm)	32.1	31.5 - 36.5
Specific Gravity:	1.280	1.26 - 1.34

**Malcolm McLaughlin, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*