



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/31/2018

CERTIFICATE OF ANALYSIS **TERGAJET**

Manufactured by: Alconox, Inc.

Sub lot No: T1G8 T1G8A, T1G8B, T1G8C, T1G8D

Sub lot numbers of Master lot No: T1M7

Manuf. Cat. no.: 2201, 2203, 2204, 2204-1, 2225, 2250

Date Of Manufacture: 7/27/2018

Expires end of: 07/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.5	11.0 - 12.5
Surface Tension:(1% solution Dyne/cm)	31.8	31.0 – 36.0
Foam (1% in DI water)	NIL	NIL

Malcolm McLaughlin, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.