



**vwr.com**  
**1.800.932.5000**

Distributed by: VWR International  
Manufacturer: Alconox, Inc.  
30 Glenn St.  
White Plains, NY 10603 USA  
Ph 914-948-4040  
cleaning@alconox.com

August 20, 2019

**CERTIFICATE OF ANALYSIS**  
**LABTONE**

Manufactured by: Alconox, Inc.  
Lot No: V1H9  
VWR Cat. No.: 21850-027, 89030-112 (4X8 lb bottle size) and 21850-080, 89030-114 (25 lb box size)  
Date Of Manufacture: 08/19/2019  
Shelf Life Expiration Date (2 years from manufacture): 08/21  
Storage Conditions: Closed upright and in a cool dry place (15 – 30 deg C or 60 – 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH: (1% solution in DI water)	9.8	9.5 – 11.5
Surface Tension (dynes/cm): (1% solution in DI water)	30.8	28.0 – 33.0

Malcolm McLaughlin