



Alconox, Inc.

The Leader in Critical-Cleaning

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 5/17/2018

## **CERTIFICATE OF ANALYSIS** **DETONOX**

**Manufactured by: Alconox, Inc.**

Sub lot No: X051618

Sub lot numbers of Master lot No: X091517

Manuf. Cat. no.: 2301, 2305, 2315, 2355

Date Of Manufacture: 5/16/2018

Expires end of: 05/20

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	Clear pale yellow liquid –nearly odorless
pH: (1% solution in DI water)	9.6	9.5 - 11.5
Surface Tension:(1% solution Dyne/cm)	32.2	29.0 - 33.0
Specific Gravity:	1.074	1.040 -1.120

**Malcolm McLaughlin, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*