

09/14/2020

## CERTIFICATE OF ANALYSIS 20% CITRIC ACID

Manufactured by: Alconox Inc.

Lot No: C091120

Manuf. Cat. no.: C201, C201-1

Date of Manufacture: 09/11/2020

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*

Shelf Life Expiration End Of: 09/2022

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Analysis</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:	OK	Clear liquid
Specific Gravity:	1.09	1.07-1.10
Citric Acid (%/wt)	20.3%	19.0-20.5%

Laura Falla, QA



Alconox Inc