



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/4/2020

CERTIFICATE OF ANALYSIS TERGAZYME

Manufactured by: Alconox, Inc.

Sub lot No: G1F0 G1F0A, G1F0B, G1F0C, G1F0D

Sub lot numbers of Master lot No: G1E9

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 6/1/2020

Expires end of: 06/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	White powder with cream colored flakes and small brown granules
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.2	9.0 - 10.0
Surface Tension:(1% solution Dyne/cm)	31.0	29.0 - 33.0

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.