



## **SOLUJET®**

**Low-Foaming Phosphate-Free Liquid Detergent**

### **Sample Trace Analysis Undiluted Concentrate**

<b>ELEMENT</b>	<b>DETECTED</b>	<b>DETECTION LIMIT</b>
<b>Chlorine</b>	NONE	Less than 100 PPM
<b>Fluorine</b>	NONE	Less than 10 PPM
<b>Bromine</b>	NONE	Less than 10 PPM
<b>Iodine</b>	NONE	Less than 10 PPM
<b>Boron</b>	NONE	Less than 10 PPM
<b>Lead</b>	NONE	Less than 10 PPM
<b>Zinc</b>	NONE	Less than 10 PPM
<b>Cadmium</b>	NONE	Less than 50 PPM
<b>Arsenic</b>	NONE	Less than 50 PPM
<b>Antimony</b>	NONE	Less than 50 PPM
<b>Bismuth</b>	NONE	Less than 10 PPM
<b>Tin</b>	NONE	Less than 10 PPM
<b>Sulfur</b>	60 PPM	Less than 10 PPM
<b>Mercury</b>	NONE	Less than 0.5 PPM
<b>Phosphorous</b>	NONE	Less than 100 PPM
<b>Potassium</b>	An ingredient in SOLUJET present at 1–5% by weight	
<b>Sodium</b>	An ingredient in SOLUJET present at 1–5% by weight	

**Note:** SOLUJET detergent is recommended for use at 1% dilutions.