



CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)

Date of Manufacture: 07/31/2020 to meet BC requirements

Manufactured by: Alconox, Inc.

Manufacturer Lot number: C1E0

BCI Lot number: C9G0

Manuf. Cat. No.: 1801B

BCI Part No: 81912

Shelf Life Expiration Date: 07/31/2022 to meet BC requirements

Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

Value	Result	Acceptance Range
Appearance:	OK	Clear pale olive colored liquid
pH: (1% solution in DI water)	2.8	2.0-3.0
Surface Tension: (1% solution Dyne/cm)	32.3	29.0 - 33.0
Specific Gravity	1.111	1.090 - 1.140
Number of Bar Code Labels Used for Labels	720	Not Applicable
Number of Bar Code Labels Used for Records	1	Not Applicable
Number of Bar Code Labels Destroyed	163	Not Applicable

Laura Falla, QA

Alconox, Inc.
The Leader in Critical Cleaning
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com

2022-07-31

REF 81912
LOT C9G0

H628819121A

+\$\$2073122C9G0AB



(01)15089590223571
(11)200731
(17)220731
(10)C9G0

2020-07-31