

CERTIFICATE OF ANALYSIS
CITRANOX
FOR BECKMAN COULTER (BCI)


Date of Manufacture: 08/31/2020 to meet BC requirements
Manufactured by: Alconox, Inc.
Manufacturer Lot number: C1G0
BCI Lot number: C9H0
Manuf. Cat. No.: 1801B
BCI Part No: 81912
Shelf Life Expiration Date: 08/31/2022 to meet BC requirements
Storage Conditions: Closed upright and in a cool dry place (15 - 30 deg C or 60 - 90 deg F).

<u>Value</u>	<u>-</u>	<u>Result</u>	<u>Acceptance Range</u>
Appearance:		OK	Clear pale olive colored liquid
pH: (1% solution in DI water)		2.9	2.0-3.0
Surface Tension: (1% solution Dyne/cm)		32.5	29.0 – 33.0
Specific Gravity		1.111	1.090 – 1.140
Number of Bar Code Labels Used for Labels		720	Not Applicable
Number of Bar Code Labels Used for Records		1	Not Applicable
Number of Bar Code Labels Destroyed		133	Not Applicable



Laura Falla, QA

REF 81912

 2022-08-31

LOT C9H0




+H628819121A



+\$\$2083122C9H0AD



(01)15099590223571
 (11)200831
 (17)220831
 (10)C9H0

 2020-08-31

101550A

The Leader in Critical Cleaning
 30 Glenn St. White Plains NY 10603 USA
 Tel 914.948.4040 • Fax 914-948-4088
 www.alconox.com • cleaning@alconox.com