



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 3/15/2021

## CERTIFICATE OF ANALYSIS ALCONOX

### Manufactured by: Alconox Inc.

Sub lot No: A3C1 A3C1A, A3C1B, A3C1C, A3C1D, A3C1E

Sub lot numbers of Master lot No: A2F0

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 3/11/2021

Expires end of: 03/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### ANALYSIS

### RESULT

### ACCEPTED RANGE

|                                                       |      |                                        |
|-------------------------------------------------------|------|----------------------------------------|
| Appearance:                                           | OK   | White powder with cream colored flakes |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete                               |
| pH:(1% solution in DI water)                          | 9.3  | 9.0 - 10.0                             |
| Surface Tension:(1% solution Dyne/cm)                 | 31.0 | 29.0 - 33.0                            |

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*