



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 12/22/2020

CERTIFICATE OF ANALYSIS DETOJET

Manufactured by: Alconox Inc.

Sub lot No: D2M0 D2M0F, D2M0G, D2M0H, D2M0J, D2M0K, D2M0L, D2M0M

Sub lot numbers of Master lot No: D1G9

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 12/15/2020

Expires end of: 12/22

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.0	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.34	1.27 - 1.34

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.