



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 6/23/2021

CERTIFICATE OF ANALYSIS LUMINOX

Manufactured by: Alconox Inc.

Sub lot No: M1F1 M1F1G, M1F1J, M1F1L, M1F1M

Sub lot numbers of Master lot No: M1C0

Manuf. Cat. no.: 1901, 1901-1, 1915, 1955

Date Of Manufacture: 6/7/2021

Expires end of: 06/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear to pale yellow liquid
pH: (as is)	6.8	6.0 - 8.0
Surface Tension:(1% solution Dyne/cm)	34.7	32.0 - 38.0
Specific Gravity:	1.02	1.00 - 1.06
Foam (1% in DI water)	NIL	NIL

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.