



Alconox, Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 7/14/2021

## CERTIFICATE OF ANALYSIS LIQUINOX

**Manufactured by: Alconox, Inc.**

Sub lot No: L1G1 L1G1F, L1G1H, L1G1J, L1G1K, L1G1L, L1G1M

Sub lot numbers of Master lot No: L1H9

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 7/13/2021

Expires end of: 07/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	Clear golden liquid
pH: (as is)	8.4	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	30.8	27.0 - 31.0
Specific Gravity:	1.08	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*